

Compressor Nebulizer

Mini-Neb

Mini-Neb is a high quality aerosol compressor with a small foot print design. The compact design and integrated handle make it easy to carry and store. With the ideal 0.5-5 micron particle, user receives an efficient and effective aerosol therapy. In addition, Mini-Neb comes with the better price package which meet the need of managed care. Portable compressor is the best type of compressors for patients and families who lead an active life and also for people with unpredictable asthma. With good medical technology base, our new Mini-Neb is your perfect solution for aerosol therapy.

Key Features:

- ◆ **Compact design**
The streamlined appearance and small size design make it easy to clean and store.
- ◆ **Consistent performance**
Reliable motor-driven piston keeps the performance consistent.
- ◆ **Comfortable and effective therapy**
Comfortable therapy is delivered with 0.5-5.0 microns particles that can go down the deeper airway easily, enhancing therapeutic effect.
- ◆ **Easy-to-carry**
A built-in handle is designed for portability even for small hands, and lets users take their aerosol treatments everywhere.



Specifications

Specification	Description
Compressor Pressure	≥40 psi
Compressor Flow	≥18 LPM
Nebulizer average rate	0.3 ml / min
Weight	2.1 kg or 4.6 lb
Size	19.6x18.5x11.6 cm or 7.7x7.3x4.6 inch
Power input	120 VAC / 60HZ, 230 VAC / 50HZ
Operating Temperature	+10°C to +40°C / +50°F to +104°F
Operating Humidity	10% - 90% RH
Storage Temperature	-15°C to +50°C / +5°F to +122°F
Storage Humidity	10% - 90% RH

* This Compressor Nebulizer is only to be used with the Salter Lab 8900

All specifications are subject to change without notice.



9, Min Sheng St. Tu-Cheng,
Taipei County, 236, Taiwan, R.O.C.
TEL: 886-2-22683100
FAX: 886-2-22686525
E-mail: ss@apexmedical.com.tw
http://www.apexmedical.com.tw