

# Compressor Nebulizer



## MiniCare

MiniCare nebulizer is a high quality aerosol compressor with contemporary, non-medical styling. It is appealing to the clinicians and the patients. Featuring tubing holder design, large storage compartment and reliable performance, it brings users convenience and a high quality aerosol therapy. With its powerful output, MiniCare delivers medications with proper size for inhalation, quickly and effectively. It is a perfect solution for those patients seeking a quality, effective treatment from a compressor nebulizer with a stylish outlook.

## Key Features:

- ◆ **Appealing outward appearance**  
The "user friendly" design of MiniCare makes children's respiratory therapy less frightening
- ◆ **Convenient built-in handle**  
The built-in handle allows the user to carry the unit easily and even move around during the treatment
- ◆ **Large storage compartment**  
Storage compartment easily accommodates nebulizer accessories and tubing after every usage
- ◆ **Hidden storage for electrical cord**  
A hidden compartment allows the electrical cord to be stored within the housing of compressor unit
- ◆ **Reliable performance**  
Equipped with reliable piston-type compressor, MiniCare has proven itself by obtaining CE/cUL/ FDA 510(k) approval



## Specifications

MiniCare	
Max. Pressure	≥ 30 psi
Max. Flow Rate	≥ 12 lpm
Operating Flow Rate @ 10 psi	≥ 5 lpm
Nebulization rate	0.2 - 0.3 ml / min*
Particle dimension	0.5 - 5 microns*
Dimensions (L x W x H)	30.0 X 16.5 X 10.0 cm / 11.8" X 6.5" X 3.9"
Weight	Approx. 1.9 kg / 4.2 lb
Electrical Requirement	120VAC / 60 Hz, 230VAC / 50 Hz
Operating Range	10°C ~ 40°C (50°F ~ 104°F) 10% ~ 90% RH non-condensing
Storage Range	-15°C ~ 50°C (5°F ~ 122°F) 10% ~ 90% RH non-condensing
EN60601-1 Classification:	Class II, Type BF, IPX0 No AP/APG Protection

\* Particle size and nebulization rate are determined by the nebulizer used with the device.

All specifications are subject to change without notice.



Private label is welcome



9, Min Sheng St. Tu-Cheng,  
Taipei County, 236, Taiwan, R.O.C.  
TEL: 886-2-22683100  
FAX: 886-2-22686525  
<http://www.apexmedical.com>