

PM210 Mattress Overlay System



PM210

PM210 is designed to offer the most trouble free and cost effective pressure relief surfaces in the health care market. The PM210 system is an alternating mattress overlay for home care use, designed to aid in the prevention and treatment of Stage I to Stage III pressure ulcers. It combines the basic but important features including static function, low pressure alarm, power failure alarm, emergency CPR and quick connector, & etc. The advanced technology of compressor operating on demand makes PM210 system more durable and prolongs its product life. PM210 offers efficient treatment and easy-to-use operation, which can provide benefits for both patients and caregivers.

Key Features:

Pump

- ◆ Easy operation pressure knob control
- ◆ Quiet, Vibration-free
- ◆ Compressor operates on demand
- ◆ Static function helps facilitating nursing care and patient repositioning
- ◆ Low Pressure Alarm with audible buzzer and visible indicator
- ◆ Power Failure Alarm to remind caregiver when power supply is interrupted
- ◆ Alarm mute to enable/disable the audible alarm

Mattress

- ◆ Low air loss feature helps to prevent heat and moisture build-up
- ◆ Emergency CPR deflation
- ◆ Transport mode by proprietary quick connector
- ◆ Vapor permeable, water resistant and fire retardant cover
- ◆ Elegant zipper-all-around cover for easy maintenance and cleaning
- ◆ Interchangeable, individual air bladder for easy maintenance



Specifications

PM210 Mattress Overlay system

Electrical	AC 120V 50/60 Hz, 0.17A AC 230V 50/60 Hz, 0.07A
Pressure	Soft to Firm setting
Air Output (Lpm)	8 liter/min
Cycle time (Mins)	Fixed, none adjustable, 8 (60Hz) / 9 (50Hz)
PUMP Dimension (L x W x H)	29.1 x 20 x 11.7 cm / 11.4" x 7.9" x 4.6"
PUMP Weight	2.2 kg / 5 lb
Mattress Model	5" x 17 Air bladders
Cover Material	Nylon/PU
Cells Material	Nylon/PU
Base Material	Nylon/PU
Mattress Dimension (L x W x H)	200 x 90 x 12.7 cm / 78.7" x 35.4" x 5" 200 x 85 x 12.7 cm / 78.7" x 33.5" x 5" 200 x 80 x 12.7 cm / 78.7" x 31.5" x 5"
Mattress Weight	5.2 kg / 11.4 lb
Weight Capacity	140kg / 300lb

All specifications are subject to change without notice.



Private label is welcome



APEX MEDICAL CORP.
9, Min Sheng St. Tu-Cheng,
Taipei County, 236, Taiwan, R.O.C.
TEL: 886-2-22683100
FAX: 886-2-22686525
<http://www.apexmedical.com>