

PM211 Mattress Replacement System



PM211

The PM211 system is a digital controlled alternating pressure relief mattress replacement for hospital bed or home use, designed to aid in the prevention and treatment of Stage I to Stage IV pressure ulcers. It comprises of many important features including pressure sensor technology, power failure alarm, static function, low pressure alarm, emergency CPR and quick connector, & etc. The pressure sensor technology of PM211 ensures a consistent, desired low-pressure support while preventing the incidents of bottoming out. The PM211 system offers variety of features that make it efficient and easy-to-use.

Key Features:

Pump

- ◆ Microprocessor controlled for simplicity, reliability and ease of use
- ◆ Quiet, Vibration-free
- ◆ Multi Therapy Modes provide patient best pressure relief and comfort - Alternating/Static/Auto-firm/Seat Inflation function
- ◆ Smart pressure detection technology increase the pump's life expectancy
- ◆ Power failure Alarm to remind caregiver when power supply is interrupted
- ◆ Low Pressure & Service Alarm to notify caregiver
- ◆ Auto-recovery static function returns the mattress back to alternating mode

Mattress

- ◆ Low air loss to provide good ventilation and keep users dry
- ◆ Emergency CPR deflation designed for quick air release
- ◆ Transport mode by proprietary quick connector
- ◆ Vapor permeable, water resistant and fire retardant cover
- ◆ Elegant zipper-all-around cover for easy maintenance and cleaning
- ◆ Interchangeable, individual air bladder for easy maintenance



Specifications

PM211 Mattress Replacement system

| | |
|--------------------------------|---|
| Electrical | AC 120V 50/60 Hz, 0.17A AC 230V 50/60 Hz, 0.07A |
| Pressure | Adjustable 10 level pressure range |
| Air Output (Lpm) | 8 liter/min |
| Cycle time (Mins) | 10, 15, 20, 25 selectable |
| PUMP Dimension (L x W x H) | 29.1 x 20 x 11.7 cm / 11.4" x 7.9" x 4.6" |
| PUMP Weight | 2.2 kg / 5 lb |
| Mattress Model | 8" x 20 Air bladders (cell-in-cell / cell-on-cell) |
| Cover Material | Nylon/PU |
| Cells Material | Nylon/PU |
| Base Material | Nylon/PU |
| Mattress Dimension (L x W x H) | 200 x 90 x 20.3 cm / 78.7" x 35.4" x 8" 200 x 85 x 20.3 cm / 78.7" x 33.5" x 8" 200 x 80 x 20.3 cm / 78.7" x 31.5" x 8" |
| Mattress Weight | 7.5 kg / 16.5 lb |
| Weight Capacity | 180kg / 400lb |

All specifications are subject to change without notice.



Private label is welcome



APEX MEDICAL CORP.
9, Min Sheng St. Tu-Cheng,
Taipei County, 236, Taiwan, R.O.C.
TEL: 886-2-22683100
FAX: 886-2-22686525
<http://www.apexmedical.com>