

Compressor Nebulizer



PenguinNeb

Today, there are over 20 million children suffered from asthma around the world. PenguinNeb's cute cartoon character is specially designed for children to minimize their negative feelings about treatment, not only can provide an effective and efficient treatment, but also can make medication fun and easy.

Key Features:

- ◆ Integrated tubing holder design holds the nebulizer kit in place easily
- ◆ The integrated handle allows easy carrying.
- ◆ The cute penguin outlook can lower children's negative feeling so that improves compliance.
- ◆ Equipped with a reliable piston-type compressor
- ◆ Optional accessories
Include a jet nebulizer, connecting tubing, mouthpiece set, pediatric and adult masks.



Easy to Carry Handle



Tubing Holder Design



Air Filter Intake

Specifications

PenguinNeb

Maximum pressure	≥ 35 psi
Maximum flow rate	≥ 13.5 lpm
Operating flow rate	5.5 lpm
Nebulization Rate	0.2 ~ 0.3 ml/min
Particle dimension	0.5 ~ 5 microns
Power Requirement	120VAC / 60Hz, 230VAC / 50Hz
Dimensions (W x D x H)	Approx. 17.5 x 12.7 x 20.3 cm / 6.9" x 5" x 8"
Weight	Approx. 2 kgs / 4.4 lbs
Operating environment	10°C ~ 40°C (50°F ~ 104°F) 10% ~ 90% RH non-condensing
Storage environment	-15°C ~ 50°C (5°F ~ 122°F) 10% ~ 90% RH non-condensing

* Particle size and nebulization rate are determined by the nebulizer used with the device.



All specifications are subject to change without notice.

Private label is welcome



APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,
Taipei County, 236, Taiwan, R.O.C.
TEL: 886-2-22683100
FAX: 886-2-22686525
<http://www.apexmedical.com>