

Continuous Positive Airway Pressure

SRT II CPAP

SRT II CPAP is designed to enhance patient's comfort and improve user compliance, and provides the ideal entry-level CPAP option for the new patient with OSA. The LCD display makes the unit convenient and easy to use for both physicians and patients, and the enhanced blower provides stable delivery of a wide range of pressures closely matched with individual's needs.

For easy hand carry, it comes with an integrated hidden handle, small size design, and an optional carrying bag. In addition, SRT II is one of the quietest CPAP devices in the market with sound level of only around 30dBA during operation. Also the total operating hours can be indicated for patient management and reimbursement. Combining many key features associated with premium CPAP devices, SRT II is sleek, quiet, and compact, making it a welcome auxiliary to the bedroom or for the traveling patient.

Key Features:

- ◆ **Easy-to-carry design**
It is compact and easy to carry with its small size and integrated hidden handle on the backside. Traveling with CPAP is no longer troublesome.
- ◆ **LCD display with back light and easy-to-use interf**
The LCD display is easy to read, and the intuitive interface is simple to use.
- ◆ **Pressure level can be adjusted based on prescription**
Operating pressure, ranging from 4 to 18cm H₂O, can be pre-set to comply with physician's prescription and closely matched with individual needs.
- ◆ **Ramp time can be set for delay delivering pressure**
Deluxe ramp capability allows a 0 - 45 minute adjustable time period to help user fall to sleep easier.
- ◆ **Compliance meter for reimbursement support**
The therapy hours can be displayed in detail at patient usage for reimbursement or therapeutic purposes.
- ◆ **Extremely quiet operation**
The delicate technical design makes it only 30dBA during operation for general use.



Specifications

SRT II CPAP

Pressure Range	4 - 18 cmH ₂ O (0.5 cmH ₂ O increment)
Ramp Time	0 - 45 mins (5 min)
Sound Level*	Under 30 dBA*
Control	START/STANDBY, RAMP, PRESSURE, METER
Indicator	LCD Display
Electrical Power	Universal power supply 100-240 VAC, 50/60 Hz Power Cord Option
Dimensions	17.1 x 25.3 x 11.0 cm / 6.7" x 9.96" x 4.3"
Weight	approx. 1.27 kg / 2.8 lb
Power Supply	100 - 240V, 50/60Hz
Operating Temperature	+5°C to +35°C (+41°F to +95°F)
Environment Humidity	15% to 95%RH non-condensing

* Measured in accordance with ISO 17510-1 at 10 cmH₂O, 1 m distance from the unit.

All specifications are subject to change without notice.

FDA **CE**
0197

Private label is welcome



9, Min Sheng St. Tu-Cheng,
Taipei County, 236, Taiwan, R.O.C.
TEL: 886-2-22683100
FAX: 886-2-22686525
<http://www.apexmedical.com>