

Continuous Positive Airway Pressure

SRT III CPAP

SRT III CPAP is the well-rounded model of the SRT CPAP series. It comes with accurate pressure, altitude compensation and detailed usage profile. Furthermore, the usage profile can be easily downloaded to PC using compatible software Easy Compliance. The software provides easy-to-read usage metr charts, daily pressure chat and statistics report of compliance data during a selected period. Printable version of records for convenient reading and filing. SRT III can store daily average pressure data and usage profile up to 365 sessions. SRT III CPAP System are user friendly, patients can enjoy the comfortable treatment with advanced features while the physicians and home care providers care providers find managing patient data very easy and efficient.

Key Features:

- ◆ **Easy-to-carry design**
Compact and lightweight, SRT III is easy to carry with integrated handle.
- ◆ **Effective treatment with altitude compensation**
Automatic leak adjustment and the manual selected altitude compensation ensure accurate output and optimal treatment effect.
- ◆ **Usage data downloadable**
With Easy Compliance Data Management Software, usage profile (compliance meter and device operating meter) can be downloaded via USB interface, allowing physician and home care provider to process further analysis.
- ◆ **Convenient charts and analysis report**
Easy Compliance provides daily pressure chart and statistics report of compliance data during a selected period, making patient management easy and efficient.
- ◆ **Data can be printed out easily**
All the records and information in the management software can be printed out for convenient reading and filing.
- ◆ **Main device can be integrated with heated humidifier (optional)**
Enhances the user comfort with moistened air delivered.
- ◆ **Optional carrying bag**



Specifications

SRT III CPAP

Pressure Range	4 - 18 cmH ₂ O (0.5 cmH ₂ O increment)
Pressure Accuracy	±0.5 cmH ₂ O
Ramp Time	0 - 45 mins (5 min)
Ramp Starting Pressure	3 -17 cmH ₂ O
Altitude Compensation	Low, Medium, High
Sound Level*	Under 30 dBA*
Dimensions	17.1 x 25.3 x 11.0 cm / 6.7" x 9.96" x 4.3"
Weight	approx. 1.29 kg / 2.8 lb
Power Supply	100 - 240V, 50/60Hz
Data Management Software	Easy Compliance
Operating Temperature	+5°C to +35°C (+41°F to +95°F)
Operating Humidity	15% to 95%RH non-condensing

* Measured in accordance with ISO 17510-1 at 10 cmH₂O, 1 m distance from the unit.

All specifications are subject to change without notice.

FDA **CE**
0197

Private label is welcome



APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,
Taipei County, 236, Taiwan, R.O.C.
TEL: 886-2-22683100
FAX: 886-2-22686525
<http://www.apexmedical.com>