

CPAP System XT Series

Continuous Positive Airway Pressure



Actual Size!

XT Series

A CPAP You Can Hold in One Hand!

XT Series CPAP System Includes:



Optional Accessories:



WiZARD-FIT™ Mask with headgear:

- ◆ Lightweight single-piece mask, made of soft medical liquid silicone, makes users feel free loading on the face
- ◆ Anatomically-contoured design with ribbed support enhances seal, comfort and fit
- ◆ 3 sizes of Full Face Mask (L, M, S) available
- ◆ 2 sizes of Nasal Mask (L, M) available



Carrying Bag:

- ◆ With soft padding provides perfect protection of CPAP device and accessories for traveling
- ◆ Inner pocket for storage of masks, tubing and power cord



Easy Compliance 2.1 Data Management Software:

- ◆ Exclusively for XT-III/XT-III Plus CPAP
- ◆ Daily pressure chart and statistics report of compliance data during a selected period makes patient management easy and efficient
- ◆ USB cable included



XT-I:

- ◆ Compact, lightweight and travel-friendly
- ◆ Entry level model for simple and economic CPAPs
- ◆ LCD display with back light and easy-to-use interface
- ◆ Operating meter directly shown on the LCD display for reimbursement support
- ◆ Ultra quiet operation



XT-III/XT-III Plus:

- ◆ Advanced model with superior features
- ◆ Enhanced pressure sensor for even better pressure accuracy
- ◆ Automatic leak adjustment and manual altitude compensation ensure accurate pressure output and optimal treatment effect
- ◆ With Data Management Software, usage profile (pressure setting, total blower hours, total therapy hours, & etc.) can be downloaded via USB interface, allowing physician and home care provider to track usage record
- ◆ Alarm clock functionality
- ◆ Premium model XT-III Plus includes AI/HI detection function



XT Heated Humidifier:

- ◆ Integrated with XT series CPAP system
- ◆ Easy to use with the rotary knob for heating level adjustment
- ◆ Water chamber with removable base plate allows ease for cleaning
- ◆ Safe & reliable with overheat protection design
- ◆ Power ON and Heater ON indicators

Features

	XT-I	XT-III/XT-III Plus
Ramp Time Function	v	v
Altitude Compensation (Manual Setting)	v	v
Leak Compensation		v
Low Pressure Auto-off		v
Data Storage	Hour Meter	Records up to 150 days usage data
Playback Software (Easy Compliance 2.1)		v
Clock/Alarm Function		v
Optional Heated Humidifier	v	v
AI/HI Detection		Only available on XT-III Plus

CPAP Technical Specifications

	XT-I	XT-III/XT-III PLUS
Pressure Range (cmH ₂ O)	4.0-18.0 (0.5 cmH ₂ O increment)	4.0-18.0 (0.5 cmH ₂ O increment)
Pressure accuracy(cmH ₂ O)	+2 -1	±0.5
Ramp Time (min)	0-45 min (5 min increments)	0-45 min (5 min increments)
Ramp Starting Pressure(cmH ₂ O)	N/A	3.0-17.0
Data Storage	Total Meter (Total Blower Hours): 9999.9 hr	Total Meter (Total Blower Hours): 9999.9 hr Compliance Meter (Total Therapy Hours): 9999.9 hr
Sound Level	under 30 dBA*	
Dimension (W x D x H)	5.8" x 5" x 3.9" (14.8 cm x 12.8 cm x 10 cm)	
Weight	1.76 lb (800 g)	
Power Supply	100V-240V, 50/60Hz, universal power supply	
Operating Range	+5°C to +35°C (+41°F to +95°F) 15% to 95%RH non-condensing	
Storage Range	-15°C to +50°C (+5°F to +122°F) 10% to 90%RH non-condensing	
Warranty	2 year	
EN60601-1 Classification:	Class I, Type BF, IPX0 IEC 60601-1-2 EMC Compatibility	

* Measured at 10 cmH₂O, 1 meter distance

Heated Humidifier Technical Specifications

	XT Heated Humidifier
Reservoir Capacity	450ml
Heater Settings:	1 to 6 from 30°C to 60°C (86°F to 140°F)
Humidity Output	10-40 mgH ₂ O/L
Pressure Drop	0.2 cmH ₂ O @ 60LPM
Dimension (W x D x H)	10" x 6.8" x 5.8" (25.5 cm x 17.2 cm x 14.6 cm)
Weight	1.37 lb (620 g)
Power Supply	120VAC, 50/60Hz (For 120V system) 220-240VAC, 50/60Hz (For 230V system)
Operating Range	+10°C to +35°C (+50°F to +95°F) 10% to 90%RH non-condensing
Storage Range	-20°C to +70°C (-5°F to +158°F) 10% to 90%RH non-condensing
Warranty	1 year
EN60601-1 Classification:	Class I, Type BF, IPX0 IEC 60601-1-2 EMC Compatibility

All specifications are subject to change without notice.



Private label is welcome