

Compressor Nebulizer

ecoNeb

ecoNeb, a small footprint unit for easy carrying to anywhere, is powerful enough to deliver a fast and effective aerosol treatment. Featuring medication cup holder design and reliable performance, it brings users ease of use and a high quality aerosol therapy.



Key Features:

- ◆ **Compact design for easy to carry**
Stylish and compact design makes it easy to clean and store.
- ◆ **Tubing holder design**
Provides an easy removable yet tight enough snug fit for the majority the tubings in the market.
- ◆ **Easy to carry**
The integrated handle allows for easy carrying.
- ◆ **Effective and efficient treatment**
Dispenses 3 c.c. medication within 10 minutes or less.
- ◆ **Optional accessories**
Includes a jet nebulizer, connecting tubing, mouthpiece, T piece, optional child and adult nasal masks, and optional carrying bag.



Specifications

ecoNeb

Max. Pressure	≥ 30 psi
Operating Flow Rate	5 lpm
Nebulization rate	0.2 - 0.3 ml / min*
Particle dimension	0.5 - 5 microns*
Dimensions (L x W x H)	Approx. 14.0 x 18.8 x 10.0 cm / 5.5" x 7.4" x 3.9"
Weight	Approx. 1.8 kg / 4.0 lb
Power Requirement	120VAC / 60 Hz, 230VAC / 50 Hz
Operating Range	10°C ~ 40°C (50°F ~ 104°F) 10% ~ 90% RH non-condensing
Storage Range	-15°C ~ 50°C (5°F ~ 122°F) 10% ~ 90% RH non-condensing

* Particle size and nebulization rate are determined by the nebulizer used with the device.



All specifications are subject to change without notice.

Private label is welcome



APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,
Taipei County, 236, Taiwan, R.O.C.
TEL: 886-2-22683100
FAX: 886-2-22686525
<http://www.apexmedical.com>