

Compressor Nebulizer



ecoPlus

The ecoPlus is a high quality, economical compressor with user-friendly design for aerosol treatment delivery.

Featuring tubing holder design and reliable performance, it brings users convenience and high quality aerosol therapy in home respiratory care.

Key Features:

- ◆ **Compact design**
Stylish and compact design makes it easy to clean and store.
- ◆ **Tubing holder design**
Provides an easy removable yet tight enough snug fit for the majority of the tubings in the market.
- ◆ **Easy to carry**
The integrated handle allows for easy carrying.
- ◆ **Effective and efficient treatment**
Dispenses 3 c.c. medication within 10 minutes or less.
- ◆ **Optional accessories**
Includes a jet nebulizer, connecting tubing, mouthpiece, T piece, optional child and adult nasal masks, and optional carrying bag.



Specifications

ecoPlus

Max. Pressure	≥ 30 psi
Max. Flow Rate	≥ 14 lpm
Operating Flow Rate @ 10 psi	5 lpm
Nebulization rate	0.2 - 0.3 ml / min*
Particle dimension	0.5 - 5 microns*
Dimensions (L x W x H)	Approx. 14.0 X 18.8 X 10.0 cm / 5.5" X 7.4" X 3.9"
Weight	Approx. 1.8 kg / 4.0 lb
Electrical Requirement	120VAC / 60Hz, 230VAC / 50Hz
Operating Range	10°C ~ 40°C (50°F ~ 104°F) 10% ~ 90% RH non-condensing
Storage Range	-15°C ~ 50°C (5°F ~ 122°F) 10% ~ 90% RH non-condensing
EN60601-1 Classification:	Class II, Type BF, IPX0 No AP/APG Protection

* Particle size and nebulization rate are determined by the nebulizer used with the device.



All specifications are subject to change without notice.

Private label is welcome



APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,
Taipei County, 236, Taiwan, R.O.C.
TEL: 886-2-22683100
FAX: 886-2-22686525
<http://www.apexmedical.com>