

# Compressor Nebulizer



## ecoStar

ecoStar, a compact unit for use at home or on the go, is a high quality, more economical compressor for aerosol treatment delivery. It's not only efficient and cost effective, but also simple to use. Featuring tubing holder design, optional carrying bag and reliable performance, it brings users convenience and a high quality aerosol therapy. ecoStar is a perfect solution for those patients seeking a quality, effective treatment from an affordable compressor with a small footprint design.

## Key Features:

- ◆ **Compact design**  
Stylish and compact design makes it easy to clean and store
- ◆ **Optimal treatment effect**  
ecoStar provides you with the efficient and effective treatment by dispensing 2-3 c.c. medication into fine particles within 15 minutes or less
- ◆ **Tubing holder design**  
Provide a convenient temporary holding for the nebulizer kit
- ◆ **Easy to carry**  
The integrated handle and lightweight design allows easy carrying to anywhere
- ◆ **Optional accessories**  
Includes a jet nebulizer, connecting tubing, mouthpiece set, optional pediatric and adult masks, and optional carrying bag
- ◆ **Reliable performance**  
Equipped with reliable piston-type compressor, ecoStar has proven itself by obtaining CE/ FDA 510(k) approval



## Specifications

ecoStar	
Max. Pressure	30 - 35 psi
Max. Flow Rate	≥ 14 lpm
Operating Flow Rate @ 10 psi	5 - 6 lpm
Nebulization rate	0.2 - 0.3 ml / min*
Particle dimension	0.5 - 5 microns*
Dimensions (L x W x H)	14.6 X 20.3 X 9.5 cm / 5.7" X 7.9" X 3.7"
Weight	Approx. 1.6 kg / 3.6 lb
Electrical Requirement	120VAC / 60 Hz, 230VAC / 50 Hz
Operating Range	10°C ~ 40°C (50°F ~ 104°F) 10% ~ 90% RH non-condensing
Storage Range	-15°C ~ 50°C (5°F ~ 122°F) 10% ~ 90% RH non-condensing
EN60601-1 Classification:	Class II, Type BF, IPX0 No AP/APG Protection

\* Particle size and nebulization rate are determined by the nebulizer used with the device.

All specifications are subject to change without notice.



Private label is welcome



## APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,  
Taipei County, 236, Taiwan, R.O.C.  
TEL: 886-2-22683100  
FAX: 886-2-22686525  
<http://www.apexmedical.com>